

# MCS<sup>®</sup>+: Multicomponent Collection System

The only true multicomponent system



**HAEMONETICS<sup>®</sup>**  
THE Blood Management Company™




# Applications

## Component donations

- Platelets
- Red blood cells
- Plasma
- All blood component combinations

## Patient treatments

- Stem cells
- Plasma exchange
- Autologous red blood cells
- Therapeutic cell removal



# Multicomponent collection

- MCS +**
- ➔ + automated Platelet Collection
  - ➔ + automated Red Cell Collection
  - ➔ + automated Plasma Collection
  - ➔ + combinations of all above
  - ➔ + therapeutic applications



# MCS<sup>®</sup>+: Multicomponent Collection System

## Key advantages

- Lightweight at 27.5 kg
- Reliable
- Good product quality
- Real multicomponent device
- Fast set-up and tear-down
- No secondary component processing
- Minimal operator involvement and time
- Quiet device



# Applications

- **Plasma** is collected to be transfused(**FFP**):
  - Substitution of coagulation factors and plasma proteins in cases of massive bleeding, transplantation, clotting disorders, burns
- **Plasma** is collected for **industrial fractionation (PPP)**:
  - Most often plasma is fractionated and sold to pharmaceutical companies as raw material for manufacture into biopharmaceuticals.
- **Platelets** are collected to be transfused:
  - Patients undergoing surgery, chemotherapy
  - Hemophilia, bleeding disorders
- **Red Blood cells** are collected to be transfused:
  - Severe hemorrhage (surgery, emergency)
  - Anemia, Thalassemia
- **Stem cells** are collected to be transfused:
  - Leukemia treatments, repopulate the whole bone marrow and grow new blood cells after chemotherapy
- **Plasma exchange** is done on a patient:
  - Removal and replacement of plasma with donor plasma or other fluids in order to treat a variety of diseases eg. Guillain-Barré syndrome or Hyperviscosity syndrome



# Plasma applications

## Plasma Collection protocols

- PPP&FFP, Platelets poor Plasma & Fresh Frozen Plasma
  - Therapeutic plasma collection with closed sets (623E, 622HS, 623HS)
  - Fractionation plasma collection with pre-connected or bundled sets (627, 792, 782, 782HS-P, 792P, 625B, 625HS, 690, 692, 699, SC690, SC692, 620, 620E)
  - Utilization of standard BMB or HS core bowl.
  - Optional Saline compensation



# Disposables


**Unbundled sets** → The customer creates his own set  
**(Commercial plasma collection):**

Needle

Donor harness

Centrifuge bowl


Collection part (**bottle or bag(s)**)



Donor harness



Bowl




Bag



# Disposables


**Pre-connected sets** → Sets are connected but not „closed“  
(PPP or FFP):





# Disposables

**Closed sets** → Sets are connected and closed: They contain bacterial filters on the connections  
**(Fresh frozen plasma collection):**



# Platelet applications

## Platelets Collection protocols

- LDP, Leukodepleted platelets (994CF-E)
  - Standard Platelets collected in plasma
  - Concurrent PPP plasma collection and Saline compensation
  - Continuous Platelet Leukocyte filtration
- LDPRBC, Leukodepleted platelets and Red Cells Collection (946FF)
  - Standard Platelets collected in plasma
  - 1 unit of RBC re-suspended in SAG-M
  - Concurrent PPP plasma collection and Saline compensation
  - Continuous Platelet Leukocyte filtration
  - Manual Red Cells leukocyte filtration



# Platelet applications



## Platelets Collection protocols

- UPP, Universal platelets protocol
  - Standard Platelets collected in plasma (997CF-E)
  - Concentrated Platelets collected in additive (999F-E)
  - Concurrent RBC re-suspended in SAG-M (949FF-E)
  - Concurrent PPP plasma collection
  - Optional Saline compensation (997CF-E)
  - Continuous Platelet Leukocyte filtration
  - Manual Red Cells leukocyte filtration (949FF-E)



# Disposables

Disposable sets for platelets collections:



# Red Cells applications

## Red Blood Cells Collection protocols

- SDR&TAE, Single Donor Red Blood Cells & Therapeutic Autologous Erythropheresis.

### SDR

- 2 units of Buffy Coat removed RBC (942)
- 2 units of Leukodepleted RBC (948F)

### TAE

- 3 units of buffy coat removed RBC (944)
  
- Red Cells re-suspended in SAG-M
- Optional Saline compensation



# Red Cells applications

## Red Blood Cells and Plasma Collection protocol


- RBCP+, Red Blood Cells & Plasma
  - 1 unit of Buffy Coat removed RBC and 2 units of plasma (941)
  - 1 unit of Leukodepleted RBC and 2 units of plasma (947F)
- Red Cells re-suspended in SAG-M
- Optional Saline compensation



# Disposables

Disposable sets for red cell collections:

**9** → MCS+ disposable set → **94XFF**  
**4** → Red cell applications → **94XFF**  
**F** → Filter for platelets → **94XFF**  
**F** → Filter for red cells → **94XFF**



# Therapeutic applications

## Therapeutic protocol


- TPE, Therapeutic Plasma Exchange protocol (980E, 981E)
  - Extract plasma from a patient
  - Substitution fluid monitored by the device
  - 2 type of bowls 125ml, 225ml
- PBSC, Peripheral Blood Stem Cells (970E, 971E)
  - Stem Cells (PBSC): collects a mononuclear cell concentrate (MNC) used for peripheral stem cell transplantations.
  - Therapeutic Platelet Reductions (TPR): reduce the platelet count of patients with abnormally high platelet counts
  - Therapeutic Leukocyte Reductions (TLR): reduce the white cell count of patients with abnormally high white cell.
  - Optional Saline compensation



# Disposables

Disposable sets for therapy:

**9** → MCS+ disposable set  
**7 or 8** → Therapy application  
**9YXE**  
**E** → International sets  
**0 or 1** → 0 for 225ml bowl  
1 for 125 ml bowl



# Questions

