

Syphilis TA EIA

An EIA for assay for the qualitative and semi-quantitative detection of IgG ,IgM & IgA antibodies to *Treponema pallidum*



CATALOGUE NUMBERS: 96 Test - DA-SYTA96
192 Test - DA-SYTA192
480 Test - DA-SYTA480

INTRODUCTION AND INTENDED USE

Syphilis is caused by the spirochaete *Treponema pallidum*, and is usually acquired by sexual contact, although the disease may be transmitted by transfusion of infected blood. Intrauterine infection also occurs. The infection is a chronic condition that typically progresses through distinct primary, secondary, tertiary, and quaternary stages of infection. These stages produce diverse clinical symptoms, typically producing initial sores known as chancres, then syphilitic rash followed by long periods of dormancy. Untreated infection may eventually result in cardiovascular problems and neurosyphilis. The organism cannot be routinely cultured in artificial media, and diagnosis of the infection usually depends on the demonstration of antibodies in the blood, which appear soon after initial infection.

Rapid Labs-Syphilis TA is for the detection of antibodies to *Treponema pallidum* in human serum and plasma, **for professional use only.**

PRINCIPLE OF THE TEST

Rapid Labs-Syphilis TA is a one step sandwich assay using recombinant antigens for enhanced sensitivity and specificity. The assay is for use with plasma and serum samples and will detect antibodies at all stages of infection. Antibodies are captured by recombinant antigens on the plate and marked by HRP recombinant conjugates for visualisation with TMB substrate. The reagents and protocol ensure ease of use and assay control.

CONTENTS

Name	Description	Colour	96 T DA- SYTA96	192 T DA- SYTA1 92	480 T DA- SYTA480
Positive Control	Human antiserum in stabilisation buffer	Red	2 mL	2 mL	2 mL
Negative control	Rabbit serum in stabilisation buffer	Yellow	3 mL	3 mL	3 mL
Plate	12 x 8 well strips coated recombinant antigen	N/A	x 1	x 2	x 5
Conjugate	HRP conjugated recombinant antigen	Orange	8 mL	16 mL	30 mL
Substrate	TMB/Peroxidase in stabilisation buffer	Pink	7 mL	14 mL	30 mL
Stop solution	0.5M Sulphuric acid	Colourless	8 mL	16 mL	30 mL
Wash buffer	20 x Concentrated	Colourless	125 mL	125 mL	250 mL

WARNINGS AND PRECAUTIONS

For *in-vitro* diagnostic use only.

Material of human origin has been tested negative by FDA Approved methods for HIV 1&2, HCV antibodies and HBsAg.

All human samples should be handled as if capable of transmitting disease and disposed of according to local guidelines.

STORAGE

Store at **2–8°C**.

Substrate is light sensitive.

LIMITATIONS OF USE

Rapid Labs-Syphilis TA may be used for neat serum and plasma.

Do not use after the stated expiry date.

Do not use substrate which has turned blue.

Controls containing sodium azide are not valid.

SAMPLES

Use fresh serum or plasma (in EDTA, sodium citrate, or heparin) samples free of microbial contamination. Samples may be stored at 2–8°C for up to 7 days prior to testing. Samples can be frozen at –20°C or lower - these should be thawed and mixed prior to testing.

ASSAY PROCEDURE

Equipment Required

Micro-pipettes capable of delivering: 50 and 300µL (+/- 10% vol)

Plate reader to read A 450nm ref 620 to 690nm and 550nm

Incubator 37°C (+/- 2°C)

Rapid Labs-SYPHILIS TA may be used in combination with automated assay equipment. Consult manufacturers for advice.

Bring all reagents and samples to room temperature before use.

Dilute wash buffer in deionised water 1/20 prior to use.

Kit controls must be run with each assay.

The kit positive control should be run in duplicate.

The kit negative control should be run in triplicate.

Sample and conjugate verification

The addition of samples, kit controls and conjugate is verified in **one step**.

The addition of conjugate to sample or kit control will cause the well to turn red, this signifies that both sample and conjugate are present. This can be confirmed by reading the wells at 550nm.

The OD will be ≥ 0.5

Substrate verification

The addition of the coloured substrate can be verified visually or by reading the wells at 550nm.

The OD will be ≥ 0.08

Protocol

1. Add 50µl neat sample or Kit Control to reaction well.
2. Add 50µl of HRP Conjugate to each reaction well, mix for 20 seconds.

Verify sample and conjugate addition as described above.

Cover the plate and incubate at **37°C for 30 minutes**.

3. Wash strips x 5 with diluted wash buffer.
Use a minimum of 300µl per wash.
A short soak of 20 seconds is recommended between washes.
Ensure excess wash is removed.

4. Add 50µl of TMB substrate to each reaction well.

Cover the plate from light and incubate at **RT for 30 minutes**.

5. Add 50µl of Stop to reaction wells.

6. Read wells at 450nm (reference filter 620–690nm)

Read within 30 minutes of addition of stop.

Assay Validation

Kit Negative:

The assay is valid if the $A_{450-620}$ value of each control reading is ≤ 0.080

If one of the values is above 0.080 then the remaining values should be used in the cut off calculation.

Kit Positive:

The value of each control reading should be ≥ 1.000

Cut-off Value

The cut off is 0.100 plus the mean of the negative values:
 $0.100 + \text{mean} (N1 + N2 + N3)$

Example: $0.100 + \text{mean} (0.020 + 0.022 + 0.021) = 0.121$

Interpretation

Negative

Samples with an OD less than the calculated cut off value are considered negative.

Positive

Samples with an OD greater than or equal to the calculated cut off value are considered positive and should be retested in duplicate.

Where both retests are below the cut off then the sample should be considered negative.

Where one of the retests is equal to or above the calculated cut off then it should be considered as a positive result and submitted for further investigation.

PERFORMANCE CHARACTERISTICS

Specificity

A study on 300 donor serum showed 100% specificity. (95% confidence limits 98.8 – 100 %)

A study on 300 donor EDTA plasma showed 100 % specificity. (95% confidence limits 98.8 – 100 %)

Clinical Sensitivity

A study on 100 positive samples showed 100 % sensitivity (95% confidence limits 96.4 – 100 %)








Analytical sensitivity

Rapid Labs-Syphilis TA has a sensitivity of 0.0015 IU/ml against the 1st IS for human syphilitic plasma IgG and IgM NIBSC code: 05/132.

BIBLIOGRAPHY

1. Larson S.A. et al: Laboratory Diagnosis and interpretation of Tests for Syphilis. Clinical Microbiology Reviews Jan.1995 Vol 8
2. Pope V. et al: A Manual of Tests for Syphilis 9TH.Edition
3. Wasley & Wong: Syphilis Serology. Oxford Medical Publications
4. Norris S.J.: Polypeptides of Treponema pallidum: Progress toward understanding their structural, functional and immunological role. Microbiological Reviews, Sept. 1993 Vol 57
5. Hagedorn H-J, Kraminer-Hagedorn A., De Bosschere K., Hulstaert F., Pottel H. and Zrein M.J. Clin. Microbiol 40 p 973-978. (2002)

KEY TO SYMBOLS

	Catalogue number
	IVD In Vitro Diagnostic Medical Device
	Manufactured by
	Temperature limitation
	Use by
	Batch code
	Consult instructions for use



Manufactured in the UK by: Rapid Labs Ltd

Unit 2 & 2A, Hall Farm Business Centre, Church Road,

Little Bentley, Colchester, Essex, CO7 8SD, U.K.

Email: info@rapidlabs.co.uk Website www.rapidlabs.co.uk