


# THE INSTRUMENTS

DAY medical's fully automated systems exist in two sizes: **DAYmate M** for medium to large numbers of samples and the **DAYmate S** as a benchtop instrument for small to medium laboratories or as a second instrument.

A streamlined user interface makes routine work smooth and simple with **minimum training time** required.

The incubation / centrifugation / reading modules can handle 1 to 18 independent tests. Incubation and centrifugation time are each standardised to 7 minutes. **Results of a full disc are available in 18 minutes independently of the profiles performed.**

Once Test Profiles are configured in our proprietary SMART software and the DAYmate consumables are loaded, minimum user intervention is required.



DAYMATE S "BENCHTOP"



DAYMATE M


→ Our range of reagents includes:

- ABO forward/reverse typing
- RhD testing
- Rh phenotyping
- K antigen determination
- Other blood group antigen typing
- Antibody screening
- Antibody identification
- DAT
- Cross matching
- Complete Quality Control kit

→ Three types of gels are available:

- DAY Neutral Gel
- DAY AHG Gel (polyspecific anti-human globulin)
- DAY AHG Gel (anti-human globulin IgG).

Reagent cells and antisera are available in dedicated cartridges sufficient for **180 tests**. Gel vials for **2400 microtubes**.



EXCLUSIVE AGENT

10, Al-Tahrir St., Sheraton Buildings, Heliopolis, Cairo, Egypt.  
 T +202 22 690 835 F +202 22 690 834 E info@curegroup.net  
[www.curegroup.net](http://www.curegroup.net)



The Innovation Setters

**DAY medical SA**

Zone Industrielle Est C5, Route de la Plaine 49, Route de la Plaine 49  
1580 Avenches, Switzerland

**T** +41 26 676 70 00 **F** +41 26 676 70 01 **E** [info@day-medical.com](mailto:info@day-medical.com)

[www.day-medical.com](http://www.day-medical.com)